ZENET

MASSAGE CAPE ZET-810

MESSAGE CAPE ZET-810

INSTRUCTIONS FOR USE

Voltage 220V

Power 9.6 W

Dear customer,

Thank you for buying one of
Models of ZENET massagers.

Please read the instruction manual and save it
for future use.

CHARACTERISTICS:

1. It will help you relax and enjoy the massage at any time that suits you.
2. Ten massage motors relax your back and hips.
3. 3-speed vibration setting.
4. Easy to use hand control with storage bag.
5. Adjustable straps are attached to your seat.
6. Soothing warming function in the lumbar spine for deeper penetration.
7. AC and DC adapters for home, car and office use.

EXPLOITATION:

1. Connect the adapter with the device to a power source.
2. Press the «POWER» button and the display lights up.
3. The «MODE» key offers 5 different programs to choose from.
4. The «HEAT» key activates the heating function and the corresponding display lights up.
5. The «INTENSITY» key to control the speed of the massage (high / medium / low).
6. Use the «TIMER» key to set the massage time (15, 30 and 60 minutes).

GUARANTEE LIABILITIES:

1. Guarantee of the device according to the operating rules - 1 year from the date of sale.
2. During the warranty period, the device must be replaced if any defects caused

by the manufacturer are found. The sales date and the stamp of the trade organization are required.

3. The guarantee does not apply to devices with:

- mechanical damage;
- defects resulting from a violation of the terms of use or incorrect actions by the owner of the device;
- deficiencies due to natural disasters;
- defects caused by the ingress of foreign bodies and liquids as a result of self-repair or self-introduction of structural changes;
- used for commercial or industrial purposes.

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